


L/S[®] Easy-Load[®] 3 Pump Heads

No tools or hardware needed

- Mount and stack heads in seconds with twist-lock mounting feature
- CRS rotor is more economical while SS offers greater durability and chemical resistance
- Mount up to four pump heads without tools or hardware when using mounting plate
- Included mounting plate adapts Easy-Load 3 pump heads to L/S drives







Pump tubing	mL per	Flow rates in mL/min		Maximum system	Catalog number‡	
sizes accepted	rev	1 to 100 rpm	6 to 600 rpm	pressure [†] psi (bar)	CRS rotor	SS rotor
Easy-Load 3 pump heads for precision pump tubing						
L/S 13	0.06	0.06 to 6	0.36 to 36	40 (2.7)		GH-77800-60
L/S 14	0.22	0.22 to 22	1.3 to 130	40 (2.7)		
L/S 16	0.8	0.8 to 80	4.8 to 480	40 (2.7)	GH-77800-50	
L/S 25	1.7	1.7 to 170	10 to 1000	35 (2.4)	<u>un-77800-30</u>	
L/S 17	2.8	2.8 to 280	17 to 1700	20 (1.4)		
L/S 18	3.8	3.8 to 380	23 to 2300	15 (1.0)		
Easy-Load 3 pump heads for high-performance precision pump tubing						
L/S 15	1.7	1.7 to 170	10 to 1000	30 (2.0)		GH-77800-62
L/S 24	2.8	2.8 to 280	17 to 1700	30 (2.0)	GH-77800-52	
L/S 35	3.8	3.8 to 380	23 to 2300	25 (1.7)	<u>un-11000-32</u>	
L/S 36	4.8	4.8 to 480	29 to 2900	20 (1.4)		
				Price		

 † Actual performance varies depending on pump tubing formulation—values shown are for firm pump tubing. [‡]CRS = cold-rolled steel SS = stainless steel

GH-77800-00 Replacement mounting kit, for mounting Easy-Load 3 pump head to Masterflex L/S pump drives. Includes mounting plate and four screws

L/S[®] Easy-Load[®] Pump Heads

Save time with fast tubing changes and reduced maintenance

- No need to remove head to change tubing
- PPS housing offers improved chemical resistance; SS rotor offers increased durability and chemical resistance
- Get four separate pumping channels by mounting up to four heads on one drive
- ATEX Zone 2 approved models rated to EEx II 3 G c II C T6
- Ready to use-include single-channel mounting hardware and 15" (38 cm) of silicone tubing







TIFIED SUPPLIER 0/516-40, -42						
log number [†]						
Black PPS ATEX-appr housing PPS hous						
SS rotor SS roto	or					

		Flour	rotoo		Catalog number [†]			
Pump tubing mL sizes per		Flow rates in mL/min		Maximum system pressure	White PPS housing		Black PPS housing	ATEX-approved PPS housing
accepted		1 to 100 rpm	6 to 600 rpm	psi (bar)	CRS rotor	SS rotor	SS rotor	SS rotor
Easy-Load pump heads for precision pump tubing								
L/S 13	0.06	0.06 to 6	0.36 to 36	40 (2.7)	GH-07516-00	GH-07516-10	GH-07518-60	GH-07516-40
L/S 14	0.21	0.21 to 21	1.3 to 130	40 (2.7)				
L/S 16	0.8	0.8 to 80	4.8 to 480	40 (2.7)				
L/S 25	1.7	1.7 to 170	10 to 1000	35 (2.4)	<u>un-0/310-00</u>			
L/S 17	2.8	2.8 to 280	17 to 1700	20 (1.4)				
L/S 18	3.8	3.8 to 380	23 to 2300	15 (1.0)				
Easy-Load pump heads for high-performance precision pump tubing								
L/S 15	1.7	1.7 to 170	10 to 1000	40 (2.7)	CH 07E1C 00	GH-07516-12	GH-07518-62	GH-07516-42
L/S 24	2.8	2.8 to 280	17 to 1700	40 (2.7)	<u>GH-07516-02</u>			<u>un-0/310-42</u>
				Price				

†PSF = polysulfone PPS = polyphenylene sulfide

CRS = cold-rolled steel SS = stainless steel



Wait! There's More at ColeParmer.com Find all the product specs you need!



Mounting Hardware

for Easy-Load Heads

Stack up to four L/S Easy-Load pump heads on one drive. Order mounting hardware based on number of mounted pump heads. Hardware is stainless steel.

Heads to be mounted	Catalog number	Price/st
One	GH-07013-04	
Two	GH-07013-05	
Three	GH-07013-08	
Four	GH-07013-09	